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
Item-no.: 143 028

Specification: Made of natural rubber, ø 25 cm, 0.6 kg, exercise instructions incl.

Important information:

This product is made of natural rubber. Slight discoloration of material are related to production process and therefore do not constitute a reason for complaint.

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A foot massage soothes not only stressed and unhappy feet, but also has essential influence on the health and well-being. Due to the spherical curvature of the Pedalo-Massage Mat, all stimulus points of the foot soles are stimulated at the same time. Standing on the Pedalo-Massage Mat for only 3 minutes each day are enough not only to release your foot musculature and consequently the deep musculature of the whole body, but also to strengthen the cross arch of foot and to prevent foot suffering. Blood circulation and purification of the legs is stimulated and the cycle is activated. This creates a pleasant feeling of warmth. In sport it is for rapid recovery from stress. The Massage Mat is ideal to carry everywhere you want and it is also rapidly applicable e. g. away on travel, after training, in the shower cubicle.

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Exercises with the Pedalo® Foot Massage Regeneration Mat

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First one, then the other foot, then both feet together. Do each of the three exercises 1 minute. Stand in all exercises with slightly bent knees, do not tense your toes and bring your ankle joints in a relaxed position.

Right foot:

Put your foot on the Massage Mat. The outside of your foot should show to the middle of the mat.

Left foot:

Repeat the exercise with your left foot.

Both feet:

Put both feet parallel on the massage mat. The insides of your foot should surround to the middle of the mat.

You can repeat the exercise, but just 1 minute per foot position! At the beginning pressure pains are normal. Always stay relaxed and let the naps act on your foot sole.

While you are sitting you can put your feet on the mat as long as it is comfortable for you. You can also sit on the mat. So your neck, shoulders, arms and hands can relax.