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Stand with both feet on the Vestimed and notice the movements of your body. Quit your tension completely and play with your balance.

Seesaw forth and back gently and relaxed by shifting weight between the forefoot and heel. Then seesaw to the left and right by weight shift to the left leg and to the right leg.

Align your body in an upright posture and hold the position.

Switch to tiptoe-position and hold the position.

Go to slight crouch through to downhill position and hold the position.

If you feel safe, perform the movements with eyes closed.

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Stand with both feet on the Vestimed and notice the movements of your body. Look at a fixed point at eye level directly in front of you.

Based on the first exercise, turn your head to your left shoulder, find your balance, then turn your head to your right shoulder, find your balance,

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then turn back to the middle and find your balance again. Alternatively with outstretched arms and a ball in one hand.

Based on the first exercise, turn your head in short quick movements several times to the left and right. Find your balance.

Based on the first exercise, move your head in circle, find your balance.

If you feel safe, perform the movements with eyes closed.

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Stand with both feet on the Vestimed and notice the movements of your body. Quit your tension completely and play with your balance.

Step into a regular rhythm on the spot. The more secure you feel, the longer you stay on one leg.

As the second exercise, additionally swing your arms.

As the third exercise change also your direction of view.

If you feel safe, perform the movements with eyes closed.

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Climb alternately slowly and evenly with one leg on the Vestimed and back down again. Gradually, increase the time spent on one leg on the Vestimed.

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As the first exercise, proceed in a series of steps over the Vestimed. While descent pay attention to do it slow.

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Stand on one foot centered on the Vestimed and keep your balance.

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Play with your free leg. Stretch it alternating times forward, sideways and backwards.

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As on the first exercise, move your foot additionally in a circle, change your direction of view.

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Go to the lateral forearm support and stabilize a straight posture.

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Stretch your free leg and arm in the extension of the body, then run elbow and knee to each other.

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Go into push-up position and keep only the upper position.

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Go into push-up position and keep only the lower position.

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Go into forearm support with your feet on the Vestimed.

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Stretch alternately one hand forward.

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As on exercise one, but forearm support on the Vestimed.

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Stand in stepping position (one leg is on the Vestimed). Perform your rearward knee toward the floor and stabilize your pelvis while doing this.

Schultergürteltraining
Training of shoulder girdle

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Stabilize your arms and shoulder girdle against the resistance of the resistance bands in different posture positions. - Arm is stretched forward / sideways

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As on exercise one, but one arm is stabilized and the second arm is leaded in different positions from bottom to top.

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Go into a slightly squat and hold the resistance bands on tension in front of your body. Perform squats.

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As the first exercise, but lead your crossed arms at the same time in front of your body from the bottom up against the resistance of the resistance bands until they have reached their full stretching.

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Lead your arms alternately in diagonal against the resistance from bottom to top.

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As the second exercise, but but when you drag with the right arm, the left leg is also splayed (one leg stand) and vice versa.

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Stand with both feet on the Vestimed (alternatively one-legged) and hold the resistance bands in front of you.

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
As on exercise one, but carry the power-bands from the front to the top in stretching.

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As exercise one, but one arm is still in front of your body (alternately move one arm).

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As you pull, stretch alternately one leg backwards. While doing this, stabilize your pelvis.

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Ganzkörperbalancetraining Whole-body balance training

One foot is secured in the attachment loop so that you can work with resistance against the resistance band. Depending to the exercise place the Vestimed, in this way, so that the resistance band is mounted in line for movement execution.

Move in the one leg stand your free leg against the resistance forward, backward, sideways from inside to outside (abduction).

As the first exercise, but you keep with the free leg opposite hand the resistance band resistance.

Schultergürtel Training Shoulder girdle training

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Stretch from the supine position the arms, which are first outstretched towards the pelvis, against the resistance of the resistance bands over your head.

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As 'shoulder gridle training' lift the pelvis additionally so that knee - pelvis - shoulder form a straight line.

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Bauch-Rumpf-Training Abdomen-trunk-training

Hold resistance bands short, lift your upper body towards your knees, hold position for a short time. Slowly quit the tension.

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Dear customer,

the Pedalo-product you have chosen is made of natural wood. Differences in color, structure, grain, small included knots or scars, prove the genuineness of our product. They are natural and therefore features of using original wood and no quality reduction.

Please note:

WARNING! The Pedalo-Vestimed is only to be used for the purposes described in this user manual.

WARNING! This device is a very demanding one! Especially when used for the first times and/or when doing exercises in standing position, a special degree in caution, attention, circumspection and a certain degree in motor abilities is assumed.

WARNING! Place the unit only on flat level ground. Protect item for damp and wet.

This device is meant for using from persons at the age of 3 up.

Maximum weight capacity approx. 120 kg.

CE Copyright by Holz-Hoerz GmbH,
printed copies also partially are not allowed.

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Like Pedalo on Facebook and always be well informed. On YouTube Pedalo-product films set you in motion and on Instagram you can see impressions from the Pedalo-theme world.