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
## Fundament unseres Körpers Fundamental base of our body

Our feet are the fundamental base of our body whose importance we mostly underestimate. The human foot is a fantastic construction of 26 bones, from many ligaments, tendons and nerves. The physician Dr. Christian Larsen discovered that nature has built our feet after the spiral principle. The back foot (heel) rotates outward, the fore foot turns oppositely inward. Therefore, the so important lengthwise and crosswise arch is arised.

If our foot arch is not intact anymore, this natural damping gets lost what results in top load on our locomotor system. A majority of the population is suffering from foot problems of which two-thirds are often caused by (wrong) movement habits and inappropriate biomechanical stress on muscles and joints, what leads to skew-, flat- or splayfeet. This inevitably means an over-stressing of foot, knee and hip joints and spine. Often this is the beginning of a chain of larger foot problems and in its further course it can lead to chronic inflammation and

increased signs of wear and tear in the joints, associated with pain. Therefore, the feet get a great significance for the whole locomotor system as basis of the straightening and body posture.

In its totality, the sole of the foot forms a sensory contact zone between body and earth. By the important perceptual signals about the nature and position of the subsoil, which it sends to the brain, the balance of the inner ear is supplemented. Also the horizontal position of the head is regulated and the sensomotoric balance of the body adjusted.

Optimal longitudinal and transverse arches of the foot are essential for a healthy foot. A healthy longitudinal and transverse arch construction provides a unique and inherently resilient construction, which allows for a statodynamic pressure distribution. In addition to the numerous foot muscles, the structures of the plantar fascia and the long soles of the foot are important factors in the development of numerous short muscles.

In addition, the interaction of the structures of the upper and lower ankle joints allows the deflection of vertical and horizontal forces acting on the foot. Further, the ankle joints ensure the mobility of the foot skeleton against the lower leg. Decisive for the function and stability of the hindfoot is the active joint stabilization and control of the heel contact during the appearance by the musculature. The round basic shape and relative small area of the heel requires the immediate balanced activation of the stabilizing lower leg and foot muscles.

Healthy feet create the basis for an optimal physical performance and maintain the necessary movement security up to the old age. With healthy feet we have fun at the movement and are able to stand, walk, run, jump, run, dance and bounce.

Begin your workout on time. We offer a high-quality and proven selection of training equipment.

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The special feature of the Pedalo 'Footboards moveable fulcrum' is the versatility. Both sides can be used and thereby enable numerous exercises for a healthy foot and a horizontally foot arch. With the flat side above, the Pedalo Footboards are comparable with the Pedalo Footboards seesaw. Thereby, strengthening and stretching of the foot as well as the longitudinal arch are the focus of attention. Turn them around, exercises are possible for the lengthwise and cross arch.

As a result, the trainees are significantly improving stability and walking safety. Tested and frequently applied in fall prevention. Nicely used for strengthening exercises in lying position when hand or forearm rests on the boards. Also ideally used as supplemental parts to balance boards or spring boards to increase level of difficulty. In combination with these additional parts, two independent movement patterns – one upon the other – have to be balanced. Stowed away quickly in every bag, this set of Pedalo-Footboards is a great companion for every training.

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#### Specification:

- 1 pair footboards in solid ash
- LxWxH: 36x14x5 cm.
- Incl. exercise instructions and anti-slip mat (48x29 cm)
- 2 kg

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For all exercises please stand relaxed, upright and with lightly bended knees. Take care about a correct foot position: Forefoot turns inwards – weight is on big toe joint. Heel turns outwards. Now adjust your leg axis by turning the lower leg inwards and the thigh outwards to the hip resulting in a straight line of ankle joint – knee joint – hip joint. Please always make sure to have a correct body posture. Stability comes out of the feet. All foot and leg muscles along the leg axis need to collude. The more safer and well controlled your stand is, the more stable your trunk, the more better your body posture and movement. While doing the exercises you will see there is a holistic effect.

**Important! Make sure that exercises are always done controlled, slowly and even fluently.**

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**Position:**

Exercise with the flat side down. Standing in lengthwise position on the Footboards. Allow stretching of the transverse arch of the foot. The foot encloses more and more the curvature of the Footboards.

**Practice:**

1. 10 sec. two-legged
2. Seesaw back and forth for 10 sec. with both legs (forefoot/heel)
3. Same as a) but each 10 sec. one-legged left/right

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
**Position:**

Exercise with the flat side down. Standing two-legged in cross position on one Footboard. Allow stretching of the transverse arch of the foot. The foot encloses more and more the curvature of the Footboards.

**Practice:**

1. Change the position of your feet several times forwards/backwards. Let each position sink in for 3 sec.
2. Lift up your legs alternately. Tread on your feet
3. Build up tension in your feet and pull your toes and heel upwards for 5 times

**Mobilisation & Stabilisation des unteren Sprunggelenks**  
Mobilization & stabilization of the lower ankle joint

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**Position:**  
Exercise with the flat side up. Standing in lengthwise position on the Footboards.

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- Practice:**
1. Standing two-legged. Hold 10 sec
  2. 5x seesawing inwards (shift your weight on the inner side of your foot) and outwards (shift your weight on the outer side of your foot)
  3. 5x seesawing two-legged in parallel left/right
  4. One-legged left/right each 10 sec

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**Position:**  
Exercise with the flat side up. Put both Footboards lengthwise one after another. Place your foot centric on both Footboards.

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- Practice:**
1. In one-leg-stand: Try to turn your foot in the metatarsus area, thus to mobilise it by turning the forefoot inwards and the heel outwards. Mobilise each foot for 15 sec.

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Mobilization & stabilization of the forefoot with fixed rear foot

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**Position:**  
The forefoot is positioned on the instable Footboard. The rear foot is positioned on the transverse Footboard.

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- Practice:**
1. Mobilization of the forefoot with conscious control
  2. Stabilize the forefoot through proprioceptive control

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Mobilization & stabilization of the rear foot with fixed forefoot

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**Position:**  
The forefoot is positioned in the middle on the transverse Footboard.

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- Practice:**
1. Turn the heel deliberately inward/outward (mobilization)
  2. Stabilize the rearfoot/heel through proprioceptive control

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Training of the shin musculature and the higher ankle joint

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**Position:**

Exercise with the flat side up. Standing two-legged in cross position on one Footboard.

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**Practice:**

1. 5x lift up your forefoot to horizontal position, keep it for 3 sec and lower again
2. 5x change from horizontal forefoot position into heel stand (lift up your forefoot), keep it for 3 sec and lower back in horizontal position. Don't put in between

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Training of the calf muscles and the higher ankle joint

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**Practice:**

1. 5x lift up your heels to a horizontal position, keep it for 3 sec and lower it again
2. 5x stand up from horizontal heel position into toe stand, keep it for 3 sec and lower in horizontal position. Don't put in between
3. Now 5x changing in an one-legged stand left/right. Hold each position for 3 sec. Don't set down your heels
4. Now 5x changing in a toe stand left/right. Hold each position for 3 sec. Don't set down your heels

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**Position:**

Exercise with the flat side down. Standing two-legged in cross position on the first Footboard.

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**Practice:**

1. Now change slowly, smooth and controlled to the second Footboard and back (left forwards, right follows, left backwards, right follows)
2. Then change the order (right forwards, left follows, right backwards, left follows)

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**Variant:**

Repeat all exercises with the flat side up.

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**Position:**

Exercise with the flat side down. Position your feet such that the toes stand in front. Now tuck up your toes forcefully against the resistance of the Footboards.

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**Performance**  
Performance

**Position:**

Exercise with the flat side up. Standing with crossed legs lengthwise on the Footboards.

**Practice:**

1. Now shift your weight to the front crossed leg. Lift your back leg and bring it in front (around the supporting leg) to return to the position on the footboard. Repeat for 3 times.
2. Same as before, but now lift the front leg and bring it backward around the supporting leg to return to the footboard.
3. Standing with both feet one after the other on the Footboards that way, that each foot is placed on one Footboard. Balance this foot position for 15 sec. Then change the foot position and repeat the exercise.

**Pedalo® Fußmassage-Regenerationsmatte**  
**Pedalo® Foot Massage Regeneration Mat**  
143028

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